

# Disposable EEG Cup Electrode Silver/Silver Chloride



**At Technomed we understand the importance of efficiency, clear signals, and ease of use.**

**Our disposable Ag/AgCl EEG Cup Electrode is designed to optimize EEG procedures and to enable you to focus on the patient and the EEG recording.**

**The electrode is designed with flexible wires and its connector is optimized for easy head box insertion and extraction.**

#### **User Benefits**

- No more cleaning and preparation of EEG electrodes
- Clear and reliable signal with optimized layer of frequency stable Ag/AgCl
- Optimized connector for easy headbox insertion and extraction
- Always new and reliable electrodes at hand
- No risk of cross contamination

## Available variants

Wire Length	Part Number
-------------	-------------

1 m	TE/C52-438
1.5 m	TE/C52-638
2 m	TE/C52-838
2.5 m	TE/C52-938

All variants are packaged in sets of 25, 5 x 5 colors: red, green, white, yellow, blue


## Specifications

### Dimensions

Diameter of cup	10 mm
Height of cup	2.6 mm

### Properties

Lead-free	✓
Latex-free	✓
BPA-free lead wire	✓
DEHP-free lead wire	✓
RoHS compliant	✓
Lead wire core	Tin plated Copper
Lead wire insulation	PVC
Connector	DIN 42802



USA
Rx ONLY



Your Technomed distributor: