

# Reusable EEG Cup Electrode Gold and Silver



**At Technomed we understand the importance of a crystal clear signal, ease of use and electrode durability.**

**Our new Reusable Au and Ag EEG Cup Electrodes are designed to give you reliable signal, and ease of use time after time.**

**Our new electrode head is designed for easy adjustment to patient contours to avoid electrode lift, and the connector is optimized for easy insertion and extraction.**

#### **User Benefits**

- Clear and reliable signal with optimized layer of Gold or Silver plating
- Elimination of electrode lift with adjustable electrode neck
- Optimal patient comfort with low profile and soft and flexible anti tangle lead wires
- Easy filling of cups with "curbed" electrode hole
- Optimized connector for easy headbox insertion and extraction
- Long electrode life with Kevlar reinforced wires

## Available variants

Materials	Wire Length	Part Number
Silver (Ag)	1 m	TE/C12-434
Silver (Ag)	1.5 m	TE/C12-634
Silver (Ag)	2 m	TE/C12-834
Silver (Ag)	2.5 m	TE/C12-934
Gold (Au)	1 m	TE/C32-434
Gold (Au)	1.5 m	TE/C32-634
Gold (Au)	2 m	TE/C32-834
Gold (Au)	2.5 m	TE/C32-934

All variants are packaged in sets of 12, 2 x 6 colors: red, green, white, yellow, blue, brown

## Specifications

### Dimensions

Diameter of cup	10 mm
Height of cup	3.1 mm

### Properties

Lead-free	✓
Latex-free	✓
BPA-free lead wire	✓
DEHP-free lead wire	✓
RoHS compliant	✓
Lead wire design	Tinsel wire
Lead wire core	Kevlar and Silver plated Copper
Lead wire insulation	TPU
Spring loaded connector	Gold plated Brass/ PVC
Connector	DIN 42802

### Cleaning and disinfection

Cleaning by washing with non-scratching tools and lukewarm water. Disinfecting with alcohol.



USA
Rx ONLY



Your Technomed distributor:

  
**Technomed**  
medical accessories

Taking you further. Step by step.

**Technomed Europe**  
Amerikalaan 71  
6199 AE Maastricht-Airport  
The Netherlands  
T +31 43 408 68 68  
info@technomed.nl  
www.technomed.nl