


# Disposable Adhesive Surface Electrode



**At Technomed we understand the importance of a crystal clear signal, ease of use, repositionability and good grounding.**

**Our new Disposable Adhesive Surface Electrode is designed with an Ag/AgCl sensor and an extra thick and flexible conductive gel. The electrodes, for recording, stimulation and ground use, are highly flexible and can be repositioned on all types of skin.**

**The electrode is designed for patient comfort and user friendliness with soft wires and an optimized connector for easy insertion and extraction.**

## **User Benefits**

- Clear and reliable signal with optimized sensor of frequency stable Ag/AgCl
- Optimal patient comfort with skin friendly soft and flexible materials
- Extra large soft and flexible gel with high repositionability
- Optimized connector for easy headbox insertion and extraction

## Available variants

Disposable 4-Disk Electrodes	Packing	Part Number
20 mm diameter disk x 4; 1 m lead wires: red, blue, green, black	10 sets	TE/K50430-001
20 mm diameter disk x 4; 2 m lead wires: red, blue, green, black	10 sets	TE/K50430-002

Disposable Ground Electrode	Packing	Part Number
40 x 50 mm surface; 1 m lead wire	20 pieces	TE/K50431-002
40 x 50 mm surface; 2 m lead wire	20 pieces	TE/K50431-003

Shelf life for unopened pouches is 24 months

## Specifications

Properties	
Lead-free	√
Latex-free	√
BPA-free lead wire	√
DEHP-free lead wire	√
RoHS compliant	√
Lead wire core	Tin plated Copper
Lead wire insulation	PVC
Connector	DIN 42802

USA
Rx ONLY



Your Technomed distributor:

